

Funky

An Unobtrusive Fingertip Health Tracker

JH3

Chuyao Feng

Hang Yang

Ze Chen

Zheng Chen


Background

- Blood carries oxygen throughout body using hemoglobin

Deoxygenated: Higher red light absorbance

Oxygenated: Lower red light absorbance

- Oxygen saturation: indicator of oxygen transport in the body




$$\text{Oxygen saturation} = \frac{C(\text{HbO}_2)}{C(\text{HbO}_2) + C(\text{Hb})} \times 100 (\%)$$

$C(\text{Hb})$ = Concentration of deoxygenated hemoglobin

$C(\text{HbO}_2)$ = Concentration of oxygenated hemoglobin



Background

- Different absorbance for deoxy Hb and oxy Hb.
- Each change in amount of red light absorbance is counted as a pulse




Design Options


- Method 1 : Transmission
- Method 2 : Reflectance
- Red and Infrared light LED are clocked to blink at different time interval
- Output signal are analyzed based on the property of photodiode



Design Pipeline




Improvable Analog Circuit Design



Transimpedance Amplifier: Convert Current Signal to Voltage Signal
Low Pass Amplifier: Cut off Frequency 10 Hz;
High Pass Amplifier: Cut off Frequency 0.5 Hz:

Block Diagram



ADC

- Ti LM331

- Sample Frequency full scale from 1KHz to 100KHz


- Low Power consumption

- 15mW at 5V max

- Adjust to 4.5mW for 10K Hz frequency.

MCU


- Ti MSP4305
- Low Power Consumption
 - 9.6mW at 3.6V (160uA/MHz)
- Dimension
 - 231.0 x 191.4 mm
- Three Channels DMA
- Two SPI & I2C bus



Estimate Cost

Item	Quantity	Price	Amount
Micro Controller	2	8.50	17.00
UR Visible Emmiter	2~3	3.15	9.45
ADC	2	3.50	7.00
Package & Design	~	50.00	50.00
Miscellaneous		50.00	50.00
Total			133.45

Software Application




Some Applications

- 1) Smart management of respiratory rehabilitation and exercise therapy.
- 2) Screening for sleep apnea syndrome.
- 3) Exercise coach.

Research stated that there would be a stratification or clustering of SpO2 % and HR as workload increased

Variables	Cluster 1	Cluster 2
SpO2 %	95 (s.d. 2.4)	84 (s.d. 4)
Heart Rate	156 (s.d. 14)	183 (s.d. 7)
Count	72	26



Q&A

Thank You